

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

### Week 4

**June 24, 2018**

Spectrophotometry Measurements at 12:00 (Stephanie)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A= 0.679
- UTEX 1% room temp. 6/13 2m A= 0.885

Spectrophotometry Measurements at 23:10 (Lin)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A= 0.890
- UTEX 1% room temp. 6/13 2m A= 0.872

Plasmid Group (Lin)

- New streaks of 1414, 1579, 2991 bacteria on spectinomycin plate

**June 25, 2018**

Spectrophotometry Measurements at 8:09 (Priya)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A= 0.870
- UTEX 1% room temp. 6/13 2m A= 0.847

Spectrophotometry Measurements at 16:00 (Natalie)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A= 0.480
- UTEX 1% room temp. 6/13 2m A= 0.493

Cell Culturing/Plating (Natalie)

- Supplementing Measurement Flasks with culture
- Added 20 mL of culture from UTEX 1% room temp. 6/13 1e into UTEX 1% room temp. 6/13 1m

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- Added 20 mL of culture from UTEX 1% room temp. 6/13 2e into UTEX 1% room temp. 6/13 2m

### Construct Group (Karthik/Matthew/Natalie/Woody/Dominika)

- G-block Resuspension of PpsbA2\_sps and PpsbA2\_cscB (followed IDT protocol)
  - Pelleted down the dry DNA for the two tubes
  - Added 25  $\mu\text{L}$  of autoclaved Milli-Q water to PpsbA2\_sps and PpsbA2\_cscB
  - Left the two tubes in the water bath at 53  $^{\circ}\text{C}$  for 15 minutes
  - Vortexed the two tubes and then centrifuged them
  - Nanodrop (first time)
    - PpsbA2\_sps= 3.7 ng/ $\mu\text{L}$
    - PpsbA2\_cscB = 4.9 ng/ $\mu\text{L}$
  - Nanodrop (second time)
    - PpsbA2\_sps= 5.1 ng/ $\mu\text{L}$
    - PpsbA2\_cscB = 8.0 ng/ $\mu\text{L}$
- Primers resuspension of 1579 PcscB start, 1579 PcscB end, 1579 Psps start, 1579 Psps end (followed IDT protocol)
  - Pelleted down the dry DNA for all four tubes
  - Resuspended 1579 PcscB start with 796  $\mu\text{L}$  nuclease free water
  - Resuspended 1579 PcscB end with 860  $\mu\text{L}$  nuclease free water
  - Resuspended 1579 Psps start with 872  $\mu\text{L}$  nuclease free water
  - Resuspended 1579 Psps end with 804  $\mu\text{L}$  nuclease free water
  - Left the four tubes in the water bath at 50  $^{\circ}\text{C}$  for 15 minutes
  - Vortexed the four tubes and centrifuged them
- PCR of PpsbA2\_sps and PpsbA2\_cscB (Short Cycle)
  - PCR Mixture 50  $\mu\text{L}$
  -

Phire Mastermix	25 $\mu\text{L}$
Forward Primer	1 $\mu\text{L}$

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification


Plasmid & Construct Design Group

Reverse Primer	1 $\mu$ L
Template DNA	2.5 $\mu$ L
H <sub>2</sub> O	20.5 $\mu$ L

- Nanodrop (Trial 1):
  - PpsbA2\_sps = 60.4 ng/ $\mu$ L
  - PpsbA2\_cscB = 40.5 ng/ $\mu$ L

### Plasmid Group (Priya/Stephanie/Manvi/Lin)

- MiniPrep x 3 (1  $\rightarrow$  NEB; 2  $\rightarrow$  Qiagen)
- Concentration of DNA = 5.4 ng/ $\mu$ L, 11.6 ng/ $\mu$ L, 58 ng/ $\mu$ L (not DNA curve)
- Attempted with Mary Lou and DNA curve  $\rightarrow$  26.7 ng/ $\mu$ L
- Split 1579, made liquid inoculation for 1414 and 2991
- [Gel extraction/gel electrophoresis for 1579](#)



- 
- Agarose gel (1%) used
- EcoRV and Sal I with cutsmart buffer
- total = 45 mL rxn, incubated 1 hour at 37  $^{\circ}$ C
- Added 9  $\mu$ L loading dye

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

**June 26, 2018**

**Construct Group (Karthik/Matthew/Natalie/Woody)**

- PCR purification of PpsbA2\_cscB (following Qiaquick protocol)
  - Didn't add sodium acetate (skipped step 1)
  - Step 7, followed the increased DNA concentration step
  - 80  $\mu$ L PB buffer
  - 16  $\mu$ L used to purify
  - 30  $\mu$ L of Milli-Q Water
  - 1  $\mu$ L used on nanodrop
  - Nanodrop of PCR purified
    - PpsbA2\_csB = 22.3 ng/ $\mu$ L
- Resuspension of PidiA\_cscB and PidiA\_sps (following IDT protocol)
  - Pelleted down the dry DNA for the two tubes
  - PidiA\_cscB  $\rightarrow$  25  $\mu$ L of autoclaved milli-Q water
  - PidiA\_sps  $\rightarrow$  25  $\mu$ L of autoclaved milli-Q water
  - Left the two tubes in the water bath at 50  $^{\circ}$ C for 15 minutes
  - Nanodrop (Trial 1):
    - PidiA\_sps = 30.3 ng/ $\mu$ L
    - PidiA\_cscB = 13.8 ng/ $\mu$ L
  - Nanodrop (Trial 2):
    - PidiA\_sps = 6.2 ng/ $\mu$ L
- PCR of PidiA\_cscB and PidiA\_sps (Short Cycle)
  - Nanodrop
    - PidiA\_cscB = 606.2 ng/ $\mu$ L(used 3  $\mu$ L)
    - PidiA\_sps= error (used 2.5  $\mu$ L)
- PCR purification of PpsbA2\_sps (following Qiaquick protocol)
  - 80  $\mu$ L of PB buffer
  - 16  $\mu$ L of PpsbA2\_sps
  - 30  $\mu$ L of autoclaved milli-Q water
  - Nanodrop (Trial 1):
    - PpsbA2\_sps = 8.2 ng/ $\mu$ L

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- Nanodrop (Trial 2):
  - PpsbA2\_sps = 7.5 ng/μL
- PCR purification of PidiA\_cscB and PidiA\_sps (following Qiaquick protocol)
  - 80 μL of PB buffer
  - 16 μL of PidiA\_cscB and PidiA\_sps
  - 30 μL of autoclaved milli-Q water
  - Nanodrop
    - PidiA\_cscB = 8.4 ng/μL
    - PidiA\_sps = 5.6 ng/μL

### Spectrophotometry Measurements at 16:08 (Lin/Natalie)

- Done at 750 nm with 1500 μL of culture
- UTEX 1% room temp. 6/13 1m A= 0.633
- UTEX 1% room temp. 6/13 2m A= 0.665

### Plasmid Group (Priya/Stephanie/Lin)

- 1579 Miniprep success
  - Nanodrop: 102.8 ng/μL
- Split 1414/2991 cultures
- Testing restriction enzyme with digest
- [Gel electrophoresis](#) → failure as EtBr was not added to solution, but found that PstI worked.

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab


Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group



- 
- 1. MW Marker (1 kb)
- 2. PstI+EcoRV
- 3. PstI
- 4. PstI + SalI
- 5. PstI + PvuII

### Cell Culturing/Plating (Natalie)

- Split the Syn. UTEX 2434 in half
- 10 mL of culture into new test tube with 10 mL of BG-11
- Supplemented the original culture with 10 mL of BG-11

### Cell Culturing/Plating (Natalie)

- Split the UTEX Collier 5/19 culture sup. With 20 mL BG-11
- 25 mL of culture into new flask labeled UTEX room temp. 6/26 (split) and supplemented it with 25 mL of BG-11
- Supplemented original culture with 25 mL of BG-11

**June 27, 2018**

### Construct Group (Karthik/Matthew/Natalie)

- Nanodrop of PCR DNA

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- PidiA\_sps = no good results
- Used 2  $\mu$ L
- PCR purify of PpsbA2\_sps
  - 80  $\mu$ L PB buffer
  - 16  $\mu$ L of PpsbA2\_sps
  - 30  $\mu$ L of EB buffer
  - Nanodrop:
    - PpsbA2\_sps = 9.9 ng/ $\mu$ L
    - Used 1  $\mu$ L
- Made .7% agarose gel
  - .4 g of agarose and 500  $\mu$ L of 100X TAE
  - Left in fridge because when pulling out comb, wells broke

### Plasmid Group (Priya/Stephanie/Sara/Lin/Manvi)

- Mini-prepped 1579
  - 102.0 ng/ $\mu$ L DNA (Had to nanodrop 2X)  $\rightarrow$  27  $\mu$ L remaining
- [Gel electrophoresis for 1579](#) with ethidium bromide staining

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab


Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group



- 
- 1. PstI + PvuII
  - Psti + PvuII, confirmation that both worked
- 2. MW Marker
- 3. PstI+EcoRV
  - Both worked
- 4. PstI + Sall
  - Both worked
- 5. PstI + PvuII
  - Did this twice because we thought this lane could have been punctured
- Gel electrophoresis w/ 1% Agarose for 1579
  - .6 g of agarose and 500  $\mu$ L of 100X TAE
  - 1. PstI + PvuII
  - 2. MW Marker
  - 3. PstI+EcoRV
  - 4. PstI + Sall
  - 5. PstI + PvuII



## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- 6 g of agarose
- Stained the gel with Diamond nucleic acid Dye that is 10,000 x concentrated
- Diluted 100X TAE into 1X TAE
- Diluted 20  $\mu$ L of 10,000X Diamond Dye with 200 mL of 1X TAE
- When transferring gel into box for shaking/staining it broke into pieces and also when transferring gel into UV box, it also broke into pieces
- Split 1579

### Spectroscopy Measurement at 17:50 (Natalie)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX Collier 5/19 Culture sup. With 20 mL BG-11 A=.697
- UTEX room temp. 6/13 (split) A= 1.009
- UTEX room temp. 6/26 (split) A= .760

### Spectroscopy Measurement at 17:56 (Natalie/Sara)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A= 0.756
- UTEX 1% room temp. 6/13 2m A= 0.751

## June 28, 2018

### Spectroscopy Measurement at 9:56 (Sara)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A=0.851
- UTEX 1% room temp. 6/13 2m A= 0.715

### Spectroscopy Measurement at 18:20 (Natalie)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A=0.933
- UTEX 1% room temp. 6/13 2m A= 0.827

Construct Group (Karthik/Natalie/Manvi/Priya/Woody)

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- Made a .7% agarose gel
  - Added .5  $\mu$ L 10,000X Diamond Nucleic Acid Dye
  - .4 g of agarose
  - 500  $\mu$ L of 100X TAE
  - Ran gel with diamond dye
    - Buffer had 18  $\mu$ L of 1X diamond dye, 3.5 mL of 100X TAE, 346.32 mL of Milli-Q water
    - Added MW marker to two wells
- Resuspension of Q3cscB and Q3sps (following IDT protocol)
  - Q3cscB  $\rightarrow$  100  $\mu$ L of autoclaved milli-Q water
  - Q3sps  $\rightarrow$  100  $\mu$ L of autoclaved milli-Q water
  - Left the two tubes in the water bath at 50  $^{\circ}$ C for 15 minutes
  - Nanodrop (used 3 $\mu$ L):
    - Q3cscB= 161.0 ng/ $\mu$ L, 92 ng/ $\mu$ L (no photos, but really smooth curves)
    - Q3cscB= 31.5 ng/ $\mu$ L
  - Nanodrop :
    - Q3sps = 292.9 ng/ $\mu$ L
- Resuspension of 1414 promoters (following IDT protocol)
  - PcpC  $\rightarrow$  50  $\mu$ L of autoclaved milli-Q water
  - PcpC560  $\rightarrow$  50  $\mu$ L of autoclaved milli-Q water
  - PidiA  $\rightarrow$  25  $\mu$ L of autoclaved milli-Q water
  - PpsbA2  $\rightarrow$  25  $\mu$ L of autoclaved milli-Q water
  - Prbc  $\rightarrow$  50  $\mu$ L of autoclaved milli-Q water
  - Left the five tubes in the water bath at 50  $^{\circ}$ C for 15 minutes
  - Nanodrop:
    - PcpC560 = 4.4 ng/ $\mu$ L, 5.4 ng/ $\mu$ L (2  $\mu$ L used)
    - PidiA = 9.5 ng/ $\mu$ L
    - PpsbA2 = 10.7 ng/ $\mu$ L
    - Prbc = 9.0 ng/ $\mu$ L, 7.7 ng/ $\mu$ L (2  $\mu$ L used)
    - PcpC = 3.4 ng/ $\mu$ L, 38.9 ng/ $\mu$ L (2  $\mu$ L used)
- Resuspension of primers (following IDT protocol)

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- 1414 promo start → 704  $\mu\text{L}$  of autoclaved milli-Q water
- 1414 promo end → 800  $\mu\text{L}$  of autoclaved milli-Q water
- 1579 cscB start → 912  $\mu\text{L}$  of autoclaved milli-Q water
- 1579 cscB end → 824  $\mu\text{L}$  of autoclaved milli-Q water
- 1579 sps start → 904  $\mu\text{L}$  of autoclaved milli-Q water
- 1579 sps end → 908  $\mu\text{L}$  of autoclaved milli-Q water
- PCR of Q3 cscB (after mixing) (CSCB cycle)
- PCR of Q3 sps (after mixing) (SPS Cycle)
- PCR of Pcp560 (EYPF Cycle)
- PCR of PpsbA2 (SHORT Cycle)

### Plasmid Group (Manvi/Priya/Stephanie)

- Miniprep with 1414 plasmid
  - Nanodrop → 203.4 ng/ $\mu\text{L}$
- Miniprep with 2991 plasmid
  - Nanodrop → 92.6 ng/ $\mu\text{L}$
- Made 1% agarose gel for staining with
  - Gel ended up breaking and having to be thrown out, might be a problem with the TAE solution

### Interlab (Lin/Matthew/Natalie)

- Finished pipetting and setting up 96 well plate for calibration 1 and calibration 2, but were unable to measure the plates because we didn't know the instrument information
- For Calibration 2, we changed tips when mixing for row E and used the same tips for the entire row for F, G and H
- Left the plate in room temperature overnight

**June 29, 2018**

### Spectrophotometer Measurement at 10:45 (Sara)

- Done at 750 nm with 1500  $\mu\text{L}$  of culture
- UTEX 1% room temp. 6/13 1m A=0.507
- UTEX 1% room temp. 6/13 2m A= 0.353

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

### Spectrophotometer Measurement at 18:28 (Lin/Lukas)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A=1.124
- UTEX 1% room temp. 6/13 2m A=0.892

### Construct Group (Woody/Matt/Natalie/Sara)

- PCR of 1414 PidiA (SHORT Cycle)
- PCR of 1414 Prbc (EYFP Cycle)
- PCR of 1414 Pcpc (EYFP Cycle)
- Resuspension of Geneblocks (following IDT protocol)
  - Q3 PisiAB\_cscB  $\rightarrow$  50  $\mu$ L of autoclaved milli-Q water
  - Q3 PisiAB\_sps  $\rightarrow$  50  $\mu$ L of autoclaved milli-Q water
  - Left the two tubes at 49  $^{\circ}$ C for 15 minutes
  - Nanodrop:
    - Q3 PisiAB\_cscB = 13.0 ng/ $\mu$ L
    - Q3 PisiAB\_sps = 11.2 ng/ $\mu$ L
- Resuspension of Geneblocks (following IDT protocol)
  - Opto EYFP  $\rightarrow$  100  $\mu$ L of autoclaved milli-Q water
  - Orig EYFP  $\rightarrow$  100  $\mu$ L of autoclaved milli-Q water
  - Left the two tubes at 50  $^{\circ}$ C for 15 minutes
  - Nanodrop:
    - Opto EYFP = 15.5 ng/ $\mu$ L
    - Orig EYFP = 8.3 ng/ $\mu$ L
- Resuspension of Q1 Primers (following IDT protocol)
  - 2991 start lone  $\rightarrow$  860  $\mu$ L of autoclaved milli-Q water
  - 2991 end lone  $\rightarrow$  912  $\mu$ L of autoclaved milli-Q water
  - Left the two tubes at 51  $^{\circ}$ C for 15 minutes
- PCR of Q3 PisiAB\_cscB (EYFP Cycle)
- PCR of Q3 PisiAB\_sps (EYFP Cycle)
- PCR of Opto EYFP (EYFP Cycle)

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- PCR of Orig EYFP (EYFP Cycle)

### Plasmid Group (Priya/Stephanie/Manvi/Sara)

- Tried making a gel with the stationary agitation method with diamond nucleic dye
  - Broke twice
- Made a gel with diamond nucleic acid dye → SUCCESS!!
  - Incubated for 20 minutes
  - Banding was smeary but indicated that ECORI and BamHi were both working

### Interlab (Lin)

- Made fluorescein stock for Calibration 3

### June 30, 2018

#### Spectrophotometer Measurement at 13:28 (Lin)

- Done at 750 nm with 1500  $\mu$ L of culture
- UTEX 1% room temp. 6/13 1m A= 1.275
- UTEX 1% room temp. 6/13 2m A= 0.984