

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

### Week 7

#### July 16, 2018

##### Interlab (Natalie/Stephanie/Lin)

- Made 180 mL of LB Media
  - 4.5 g of LB powder
  - 180 mL of milli-Q water
  - 180  $\mu$ L of chloramphenicol
- Put around 10 mL per falcon tube
- Transferred 2 colonies from each plate into liquid media
- Incubated the colonies overnight at 37 °C and 220 rpm

##### Constructs Group (Karthik/Matthew)

- Made a gel with the new promega, the bands showed up
- Tried to do a gel purification, but the gel slice wasn't melting

##### Plasmid Group (Stephanie/Sara)

- Minipreped, RE Digest, DNA purification of pAM2991 and pAM1414

#### July 17, 2018

##### Interlab (Natalie/Sara/Lin)

- Made 270 mL of LB Media
  - 6.75 g of LB powder
  - 270 mL of milli-Q water
  - 270  $\mu$ L of chloramphenicol
- Diluted the overnight cultures from 7/16
  - .5 mL of culture
  - 4.5 mL of LB media
  - Took absorbance readings and saw that the E. Coli were not growing enough
    - which means need to redo Day 2 incubation
    - Mary Lou said that we tightened the caps too much and suffocated the cells

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

-

	Colony 1	Colony 2
Positive control	.044	.044
Negative Control	.041	.043
Test Device 1	.043	.044
Test Device 2	.045	.046
Test Device 3	.045	.045
Test Device 4	.045	.043
Test Device 5	.043	.044
Test Device 6	.050	.043
LB and CAM=	.039	

-

- Redid day 2
  - Added 10 mL of LB and chloramphenicol to 16 tubes
  - Transferred 2 colonies from each plates into liquid media
  - Incubated the colonies overnight at 37 °C and 220 rpm

### Plasmid Group (Stephanie/Sara)

- Miniprep, RE Digest, DNA purification of pAM2991 and pAM1414

### July 18, 2018

#### Interlab (Natalie/Lin)

- Made 270 mL of LB Media
  - 6.75 g of LB powder
  - 270 mL of milli-Q water
- Measured the absorbance of the cultures

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

### - Colony 1

Sample	Abs600 Reading	Preloading culture (μL)	Volume of preloading media (μL)
Positive control	.153	2162	9838
Negative control	.164	1967	10033
Test device 1	.114	3333	8667
Test device 2	.165	1951	10049
Test device 3	.155	2124	9876
Test device 4	.129	2759	9241
Test device 5	.108	3636	8364
Test device 6	.159	2051	9949
Media + CAM	.042		

### - Colony 2

Sample	Abs600 Reading	Preloading culture (μL)	Volume of preloading media (μL)
Positive control	.149	2182	9818
Negative control	.150	2162	9838
Test device 1	.109	3429	8571
Test device 2	.132	2581	9419
Test device 3	.147	2222	9778
Test device 4	.112	3288	8712

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group


Experimental Verification

Plasmid & Construct Design Group

Test device 5	.096	4211	7789
Test device 6	.156	2051	9949
Media + CAM	.039		

### Constructs Group (Karthik/Matthew)

- Made a gel of DNA ladder with promega agarose, soaked in fresh diamond dye
  - Purified the gel slices using Dr. Gergen's 16000 x g centrifuge
  - Very bad DNA curves on the Nanodrop
- Made a new gel of Q1 lone and combo cscB (Phire PCR) and Q3 cscB and Q3 sps (Q5 PCR) with promega agarose, soaked in fresh diamond dye (1 or 2 bands indicates that gel purification is not required and only PCR purification is necessary)



- - Lane 1: Promega 1 kb ladder
  - Lane 2: Q3 sps PCR product from Q5 polymerase
  - Lane 3: Q3 cscB PCR product from Q5 polymerase
  - Lane 4: Q1 combo cscB PCR product from Phire polymerase
  - Lane 5: Q1 lone cscB PCR product from Phire polymerase
- Conclusions: Phire > Q5? (But for temps during PCR annealing temp for Phire set below calculated and annealing temp for Q5 set at calculated, so could be a factor)

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

**July 19, 2018**

Interlab (Natalie/Manvi/Sara/Dominika/Lin)

- Made 16 agar plates with 130 mL of LB and agar
  - Added 3.25g of LB media and 1.56g of Agar
  - 130  $\mu$ L of chloramphenicol
  - 130 mL of milli-Q water
  - Put in around 10 mL per agar plate
- Made 170 mL of LB Media
  - 4.25 g of LB powder
  - 170 mL of milli-Q water
  - 170  $\mu$ L of chloramphenicol
  - Put 10 mL of media into 16 tubes
  - Transferred 2 colonies from each agar plate into liquid media
  - Incubated the colonies overnight at 37 °C and 220 rpm
- 2 tubes spilled in the incubator, so redid them with 20 mL of LB media
  - .5 g of LB powder
  - 20 mL of milli-q water
  - 20  $\mu$ L of chloramphenicol
  - Put 10 mL of media into 2 tubes
  - Transferred 2 colonies from each plates into liquid media
  - Incubated the colonies overnight at 37 °C and 220 rpm

Plasmid Group (Stephanie/Sara)

- Miniprep, RE Digest, DNA purification of pAM2991 (realized stock on Monday/Tuesday did not have the correct enzymes)

Constructs Group (Karthik/Matthew/Woody)

- PCR adjustments to Q3 cscB and Q3 sps:
  - 1. Raise lid temp to 99, raise annealing temp to calculated, number of cycles at 30
  - 2. Touchdown: same fixes but also start annealing at 10 above and decrease by 1 each cycle till calculated annealing temp reached

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- Run on gel and soaked in Diamond dye: No discernable difference b/w 1 and 2, same non-specificity as before (multiple bands)
- Second PCR adjustments to same constructs:
  - 3. Same as 2 (with touchdown) but use 1.5  $\mu\text{L}$  template instead of 2.5, chill on ice and put in thermocycler when machine at  $>95$  during preheating
  - No gel run
- Gel purification of NEB ladder using EtBr gel and both kits
  - For both kits, used 500 and 1500 bp bands
  - Used faster centrifuge at 16,000g
  - None of the purifications worked well according to nanodrop curves

### Cell Culture/Plating (Elon)

Unless otherwise stated, for all measurements, assume an uncertainty of  $\pm 10 \mu\text{moles/s/m}^2$

- PAR Measurements Room-light:
  - 2 lights on:  $0.22\text{-}0.24 * 10^{16}$  photons/sec/cm<sup>2</sup>:  $36.52 \mu\text{moles/s/m}^2$
  - 1 light on:  $0.16\text{-}0.18 * 10^{16}$  photons/sec/cm<sup>2</sup>:  $39.84 \mu\text{moles/s/m}^2$
- Shaker Tray:
  - $0.56 * 10^{16}$ ,  $0.76 * 10^{16}$ ,  $0.62 * 10^{16}$  photons/sec/cm<sup>2</sup>
    - Range:  $92.96 - 126.16 \mu\text{moles/s/m}^2$
- Effect of papers:
  - No paper:  $0.56 * 10^{16}$  photons/sec/cm<sup>2</sup>:  $92.96 \mu\text{moles/s/m}^2$
  - Printer Paper:  $0.18 * 10^{16}$  photons/sec/cm<sup>2</sup>:  $29.88 \mu\text{moles/s/m}^2$ 
    - Dimming effect:  $63.08 \mu\text{moles/s/m}^2$
  - Parchment/Wax Paper:  $0.48 * 10^{16}$  photons/sec/cm<sup>2</sup>:  $79.68 \mu\text{moles/s/m}^2$ 
    - Dimming effect:  $13.28 \mu\text{moles/s/m}^2$
- Distance vs Light intensity (4 fluorescent lights):

Distance (inch)	Distance (cm)	Light Intensity (center)	Light Intensity (edge)	Moles of photons (center)	Moles of photons (edge)
1	2.54	2.30E+16	1.80E+16	381.80	298.80
2	5.08	1.90E+16	1.60E+16	315.40	265.60

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

3	7.62	1.70E+16	1.45E+16	282.20	240.70
4	10.16	1.50E+16	1.20E+16	249.00	199.20
5	12.7	1.40E+16	1.00E+16	232.40	166.00
6	15.24	1.20E+16	9.00E+15	199.20	149.40
7	17.78	1.00E+16	8.00E+15	166.00	132.80
8	20.32	9.60E+15	7.60E+15	159.36	126.16
9	22.86	8.50E+15	7.00E+15	141.10	116.20
10	25.4	7.90E+15	6.80E+15	131.14	112.88
11	27.94	7.20E+15	6.30E+15	119.52	104.58
12	30.48	6.70E+15	5.70E+15	111.22	94.62
13	33.02	6.20E+15	5.60E+15	102.92	92.96
14	35.56	5.60E+15	5.20E+15	92.96	86.32
15	38.1	5.30E+15	4.80E+15	87.98	79.68
16	40.64	4.60E+15	4.60E+15	76.36	76.36

- Troubleshooting in CO2 incubator: (note: by red lights I mean warm white lights, and by blue lights I mean cool blue lights)
  - Tray 3 from bottom:
    - Red lights full power:  $0.12\text{-}0.14 \times 10^{17}$  photons/sec/cm<sup>2</sup>: 199.20 - 232.40  $\mu\text{moles/s/m}^2$
    - Red & Blue full:  $0.34 \times 10^{17}$  photons/sec/cm<sup>2</sup>: 564.40  $\mu\text{moles/s/m}^2$ . Can range from 0.32 - 0.36 at edges: 531.20  $\mu\text{moles/s/m}^2$  - 597.60  $\mu\text{moles/s/m}^2$
    - Blue full:  $0.21 \times 10^{17}$  photons/sec/cm<sup>2</sup>: 348.60  $\mu\text{moles/s/m}^2$ . Can range from 0.2-0.23 at edges: 332.00  $\mu\text{moles/s/m}^2$  - 381.80  $\mu\text{moles/s/m}^2$
  - Tray 3 from bottom, but with another tray on the top slot, blocking a decent amount of light:
    - Full power on both:  $0.1 \times 10^{17}$  photons/sec/cm<sup>2</sup>: 166.00  $\mu\text{moles/s/m}^2$ . ranges from 0.9 to 1.1: 149.40  $\mu\text{moles/s/m}^2$  - 182.60  $\mu\text{moles/s/m}^2$

## COLOR CODING KEY

Spectrophotometry Measurements

Construct Group

Plasmid Group

Interlab

Cell Culture/Plating

Biobrick Group

Cyanobacteria Transformation Group

Experimental Verification

Plasmid & Construct Design Group

- Same as before, but PWM for blue lights at 155 and PWM for red at 255 (full power):
  - Tray 3:  $0.8 \pm 0.1 \times 10^{16}$  photons/sec/cm<sup>2</sup>:  $132.80 \pm 16.6$   $\mu$ moles/s/m<sup>2</sup>
  - Top tray:  $0.8 \pm 0.2 \times 10^{17}$  photons/sec/cm<sup>2</sup>:  $1,328.00 \pm 332$   $\mu$ moles/s/m<sup>2</sup>

**July 20, 2018**

Constructs Group (Karthik/Matthew)

- 3rd PCR from previous day run on gel
  - Bands very distinct, generally pure
  - Conclusion: using less template and adjusting procedure worked
- Another PCR run and examined on gel
  - Adjustments: even less template (1  $\mu$ L template)
  - Gel result: bands too faint
  - Conclusion: 1  $\mu$ L (10 ng) template not enough for PCR

Interlab (Natalie/Lin/Stephanie/Matthew/Karthik)

- Made 150 mL of LB Media
  - 3.75 g of LB
- Took overnight cultures out of the incubator, diluted them and completed day 2 for interlab
- Made 25 agar plates with 210 mL of LB and agar
  - Added 5.25 g of LB media and 2.52 g of Agar
  - 210  $\mu$ L of chloramphenicol
  - 210 mL of milli-Q water
  - Put in around 10 mL per agar plate

Plasmid Group (Stephanie/Sara)

- Miniprep, RE Digest, DNA purification of pAM1579